


# Drying Solvent Extracts using DryDisk® Membrane

William R. Jones and Zoe A. Grosser

## Key Words

Drying, Solvent Extracts, solid phase extraction, SPE, Water Extraction, Membrane, DryDisk, DryDisk-R



## Introduction


When analytes are extracted from aqueous samples with solvent either in a liquid-liquid (LLE) format or using solid phase extraction (SPE) it is likely that a small amount of water will be carried into the extract. This water should be removed before the analytical step to ensure that there is no back extraction of analytes into the water and the water is not available to cause damage to the chromatography system.

Removing water with a membrane rather than the older technique of passing the solvent through a column of sodium sulfate ( $\text{Na}_2\text{SO}_4$ ) brings several advantages. The most important analytically is that the membrane will not adsorb analytes or contaminate the extract with matrix or other potential interferences.

It is important that the membrane act effectively to keep the water on one side of the membrane and pass the solvent through.

Several experiments were performed to ensure proper operation of the DryDisk and DryDisk-R membranes. The first experiment was designed to evaluate the capacity of the membrane to separate water as a larger percentage of intermediate polar solvent, one that is both soluble in the non-polar extraction solvent and the water, was added to the solvent mixture. Water breakthrough was evaluated by looking for droplets of colored water pulled through the membrane (Figure 1).

Conditions: 30 mL of a Dichloromethane (DCM):acetone solvent mix, with 5 mL of



*Figure 1: Water (with Blue Color) and Solvent, for Separation*

Table 1: Water Breakthrough with Varying Water-Soluble Solvent Ratios

| DCM (%) | Acetone (%) | Water (mL) | Breakthrough |
|---------|-------------|------------|--------------|
| 100     | 0           | 5          | No           |
| 90      | 10          | 5          | No           |
| 80      | 20          | 5          | No           |
| 70      | 30          | 5          | No           |
| 60      | 40          | 5          | No           |
| 50      | 50          | 5          | No           |
| 45      | 55          | 5          | Yes          |

deionized water is pulled through the hydrophobic membrane. Table 1 shows the results.

This is an example of a range of solvent mixtures that may be encountered in developing an extraction process and the DryDisk-R performed well. These values will vary depending upon the exact conditions, so we recommend using less than the maximum acetone to avoid breakthrough.

A second experiment was performed to see if both types of the DryDisk gave satisfactory recoveries of a full list of analytes that might be extracted into a solvent such as DCM. The conditions of the experiment are as follows:

#### Experimental Conditions

|                |   |
|----------------|---|
| Solvent:       | 25% ethylacetate: 75% DCM. Total volume 30 mL   |
| DryVap Setting | Vacuum -10 in Hg, Nitrogen Pressure set at 27 psi. , Heat power 5, 5 µg spike of 525.2 analytes |

The recoveries , shown in Table 2, are excellent and don't significantly differ from one type of DryDisk to the other.

Table 2: Recoveries using DryDisk Membranes for Drying

| Analyte                   | #  | DryDisk R | DryDisk | Analyte                    | #   | DryDisk R | DryDisk |
|---------------------------|----|-----------|---------|----------------------------|-----|-----------|---------|
| Isophorone                | 1  | 97.0      | 100.0   | Cyanazine                  | 53  | 104.0     | 101.4   |
| 2-Nitro-m-xylene          | 2  | 101.0     | 105.2   | Metolachlor                | 54  | 101.6     | 101.4   |
| Naphthalene               | 3  | 100.2     | 104.6   | Chlorpyrifos               | 55  | 99.0      | 100.2   |
| Dichlorvos                | 4  | 99.6      | 99.8    | Aldrin                     | 56  | 100.4     | 102.0   |
| Hexachlorocyclopentadiene | 5  | 86.4      | 85.8    | Triadimefon                | 57  | 104.8     | 102.4   |
| EPTC                      | 6  | 99.0      | 101.4   | Dacthal                    | 58  | 102.8     | 104.2   |
| Mevinphos                 | 7  | 97.8      | 98.0    | MGK-264-A                  | 59  | 103.6     | 100.6   |
| Butylate                  | 8  | 94.4      | 103.4   | Diphenamid                 | 60  | 106.8     | 102.8   |
| Vermolate                 | 9  | 101.6     | 101.4   | MGK-264-B                  | 61  | 107.6     | 101.8   |
| Dimethyl phthalate        | 10 | 93.8      | 103.0   | Heptachlor epoxide B       | 62  | 101.6     | 101.2   |
| Pebulate                  | 11 | 107.6     | 105.2   | Heptachlor epoxide A       | 63  | 97.6      | 101.2   |
| Etridiazole               | 12 | 98.8      | 98.2    | Fluoranthene               | 64  | 107.0     | 104.0   |
| 2,6-Dinitrotoluene        | 13 | 103.2     | 102.4   | g-Chlordane (trans)        | 65  | 93.2      | 101.6   |
| Acenaphthylene            | 14 | 98.8      | 103.4   | Stirofos                   | 66  | 98.6      | 101.8   |
| Chloroneb                 | 15 | 102.2     | 104.0   | Disulfoton sulfone         | 67  | 97.2      | 98.8    |
| Tebuthiuron               | 16 | 100.6     | 90.4    | Butaclor                   | 68  | 98.0      | 101.2   |
| 2,4-Dinitrotoluene        | 17 | 100.4     | 100.6   | a-Chlordane (cis)          | 69  | 100.4     | 103.6   |
| Molinate                  | 18 | 101.8     | 104.2   | Endosulfan I               | 70  | 97.4      | 105.6   |
| Diethyl phthalate         | 19 | 118.6     | 106.2   | Pyrene-d10                 | 71  | 98.8      | 103.2   |
| Fluorene                  | 20 | 113.0     | 103.4   | Pyrene                     | 72  | 97.4      | 106.2   |
| Propachlor                | 21 | 112.6     | 103.6   | Napropamide                | 73  | 93.4      | 99.6    |
| Ethoprop                  | 22 | 111.8     | 103.0   | trans-Nonachlor            | 74  | 90.8      | 101.4   |
| Cycloate                  | 23 | 105.0     | 102.8   | 4,4'-DDE                   | 75  | 98.4      | 109.0   |
| Chlorpropham              | 24 | 100.4     | 100.0   | Dieldrin                   | 76  | 95.6      | 105.2   |
| Trifluralin               | 25 | 113.2     | 93.0    | Tricyclazole               | 77  | 101.0     | 86.6    |
| a-BHC                     | 26 | 100.6     | 103.0   | Terphenyl-d14              | 78  | 100.6     | 104.1   |
| Atraton                   | 27 | 99.8      | 94.0    | Endrin                     | 79  | 103.0     | 108.8   |
| Hexachlorobenzene         | 28 | 94.8      | 100.6   | Chlorobenzilate            | 80  | 95.4      | 98.2    |
| Prometon                  | 29 | 104.2     | 102.8   | Endosulfan II              | 81  | 101.8     | 110.4   |
| Lindane (g-BHC)           | 30 | 104.0     | 112.0   | 4,4'-DDD                   | 82  | 93.2      | 98.0    |
| Simazine                  | 31 | 102.8     | 102.6   | Endrin Aldehyde            | 83  | 96.8      | 102.4   |
| Atrazine                  | 32 | 97.4      | 100.8   | Butyl benzyl phthalate     | 84  | 100.4     | 104.8   |
| Propazine                 | 33 | 100.4     | 102.6   | Norflurazon                | 85  | 104.2     | 103.4   |
| b-BHC                     | 34 | 104.4     | 112.0   | 4,4-DDT                    | 86  | 93.4      | 98.4    |
| Pentachlorophenol         | 35 | 105.6     | 98.8    | Endosulfan Sulfate         | 87  | 100.4     | 106.8   |
| Pronamide                 | 36 | 100.6     | 101.4   | Bis(2-ethylhexyl)adipate   | 88  | 100.4     | 102.0   |
| Diazinon                  | 37 | 99.8      | 101.2   | Hexazinone                 | 89  | 98.8      | 103.2   |
| d-BHC                     | 38 | 109.4     | 111.8   | Triphenylphosphate         | 90  | 101.0     | 101.0   |
| Phenanthrene              | 39 | 100.0     | 102.8   | Endrin Ketone              | 91  | 96.2      | 102.6   |
| Methyl paraoxon           | 40 | 99.8      | 111.0   | Methoxychlor               | 92  | 99.2      | 97.4    |
| Anthracene                | 41 | 107.6     | 104.8   | Benz(a)anthracene          | 93  | 107.4     | 102.8   |
| Terbacil                  | 42 | 104.8     | 103.2   | Chrysene                   | 94  | 107.4     | 108.0   |
| Chlorothalonil            | 43 | 101.0     | 100.4   | Bis(2-ethylhexyl)phthalate | 95  | 104.8     | 103.0   |
| Metribuzin                | 44 | 100.0     | 93.0    | Fenarimol                  | 96  | 99.4      | 95.4    |
| Simetryn                  | 45 | 100.8     | 89.4    | cis-Permethrin             | 97  | 100.8     | 105.4   |
| Heptachlor                | 46 | 98.2      | 100.2   | trans-Permethrin           | 98  | 103.6     | 97.4    |
| Ametryn                   | 47 | 97.2      | 88.6    | Di-n-octyl phthalate       | 99  | 103.6     | 100.6   |
| Alachlor                  | 48 | 100.6     | 97.2    | Benzo(b)fluoranthene       | 100 | 114.0     | 102.6   |
| Prometryn                 | 49 | 100.0     | 94.4    | Benzo(k)fluoranthene       | 101 | 111.2     | 102.0   |
| Terbutryn                 | 50 | 99.4      | 88.6    | Benzo(a)pyrene             | 102 | 108.2     | 100.2   |
| Di-n-butyl phthalate      | 51 | 102.4     | 105.6   | Fluridone                  | 103 | 97.8      | 94.6    |
| Bromacil                  | 52 | 86.0      | 99.8    | Perylene-d12               | 104 | 98.6      | 97.0    |
|                           |    |           |         | Indeno(1,2,3-cd)pyrene     | 105 | 104.4     | 108.0   |
|                           |    |           |         | Dibenz(ah)anthracene       | 106 | 109.0     | 99.2    |
|                           |    |           |         | Benzo(ghi)perylene         | 107 | 109.4     | 103.8   |

An easier way to visualize this data to see the graphic representation of the compounds listed in Table 2. The graph below shows the recoveries in Figure 2.




Figure 2. Recoveries of a large suite of compounds after drying with DryDisk

One concern when using a drying material is if the material will introduce contamination into the extract. A test was done to see if analytes of interest spiked into an extract at low concentrations, in fact yielded more than 100% recovery.<sup>1</sup>

## Procedure

### Complete Assembly Test

1. Prepare two 1000-mL reagent water aliquots
2. Adjust the samples' pH to be less than 2
3. Spike the complete analyte mix into the samples
4. Extract the samples with conventional LLE
5. Dry one sample with sodium sulfate
6. Dry the other sample with DryDisk
7. Analyze both samples with GC/MS and compare results

The samples were spiked with a full semivolatile (method 8270) mix but, as most contamination typically appears in the form of phthalates, the data in Table 4 was condensed to focus on these compounds. As seen from the data, the recoveries for the various phthalate compounds closely match those from the sodium sulfate run. This indicates that the DryDisk membrane and holder assembly introduced no additional phthalates.

Table 4: Contamination Evaluation, Analyte Percentage Recoveries

| Analytes                    | DryDisk    |            | Na <sub>2</sub> SO <sub>4</sub> |            |
|-----------------------------|------------|------------|---------------------------------|------------|
|                             | Conc (ugL) | % Recovery | Conc (ug/L)                     | % Recovery |
| Phenol                      | 6.05       | 37.8       | 6.10                            | 38.1       |
| Napththalene                | 12.03      | 75.2       | 12.74                           | 79.6       |
| 2,6-Dichlorophenol          | 12.23      | 76.4       | 12.97                           | 81.1       |
| Dimethylphthalate           | 12.80      | 80.0       | 17.20                           | 107.5      |
| 2,4-Dinitrophenol           | 8.43       | 52.7       | 6.39                            | 39.9       |
| Pentachlorobenzene          | 12.64      | 79.0       | 13.22                           | 82.6       |
| 4-Nitrophenol               | 14.31      | 89.4       | 13.92                           | 87.0       |
| 2,3,4,6-Tetrachlorophenol   | 13.86      | 86.6       | 13.81                           | 86.3       |
| Diethylphthalate            | 13.50      | 84.4       | 13.31                           | 78.6       |
| Pentachlorophenol           | 13.80      | 86.3       | 12.57                           | 78.6       |
| Methylparathion             | 16.20      | 101.3      | 15.61                           | 97.6       |
| Heptachlor                  | 13.92      | 87.0       | 14.39                           | 89.9       |
| Di-n-butylphthalate         | 16.26      | 95.4       | 16.35                           | 102.2      |
| Aldrin                      | 14.30      | 89.4       | 14.65                           | 91.6       |
| Bis(2-ethylhexyl) adipate   | 14.93      | 93.3       | 15.66                           | 97.9       |
| Butyl benzyl phthalate      | 15.13      | 94.6       | 15.81                           | 98.8       |
| Chrysene                    | 14.86      | 92.9       | 14.94                           | 93.4       |
| 3, 3'-Dichlorobenzidine     | 13.70      | 85.6       | 13.84                           | 86.5       |
| Bis(2-ethylhexyl) phthalate | 16.55      | 103.4      | 15.34                           | 95.9       |
| Di-n-octyl-phthalate (CCC)  | 14.87      | 92.9       | 15.06                           | 94.1       |

## Conclusions

The tests done here show the DryDisk membranes perform very well in all the aspects considered. They remove water, even when a significant percentage of a water soluble solvent is mixed with the nonsoluble solvent used for extraction. In addition, the analytes show excellent recovery with both the DryDisk and DryDisk-R. The DryDisk was further compared to sodium sulfate to see if recoveries could be distinguished and base compounds pass through the DryDisk, while they show significant loss on sodium sulfate.

DryDisk membranes are a useful way to quickly remove water from organic extracts and do not add contamination or retain analytes. In addition, because the separation is physical and not chemical the amount of water that can be removed is unlimited.

## Reference

1. Susan Petitti, A Study of DryDisk Background Contamination, AN051-091214.